




NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

WIDTH EXTENSION KIT / IZYCO			
Scale: 1:8	Drafted	Verified	Confirmed
Date	13.11.2019	27.05.2025	27.05.2025
Initials	DR	KK	AZ

XXX
Plan No. E00106144/A1.1


- Use only with metal side rails (EP00256 & EP00257).
Do not use with wooden bars (EP00280) and SSR (EP00493)

Material			
No.	Plan	Article	Quantity
1#	YP08882	2629303	4
2	YP07849	26005...	4
3	P00106145	29087...	1
4	P00106146	29088...	1
5	P00106147	29089...	1
6	P00106148	29090...	1
7#	XP01795	2993539	4
8#	XP02992	GP22400	4
9#	XP01899	GP12700	8
10#	XP02466	GP21300	8
11#	M6x25	1001494	8
12#	M6x16	1008459	16(#)/8
13	Velcro 25mm	1990663	0,92m
14	Cap 50x20	1998691	8
15#	Spanring DIN471 15,7	2993606	4
16	XP01912	2993607	4
17	XP03012	2996629	4
18	XP03013	2996630	4
19	Ø5x40 DIN7	2994932	4
20#	Rubber cap Ø10	1900072	8

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

WIDTH EXTENSION KIT / IZYCO				Plan No.
				E00106144/A1.1
Scale:	1:8	Drafted	Verified	Confirmed
Date	13.11.2019	27.05.2025	27.05.2025	
Initials	DR	KK	AZ	

