



37
Detail A
Scale: 1:2

EXPLODED VIEW 276				Plan No.
				1/2 E000994/H
Scale: 1:8	Drafted	Verified	Confirmed	
Date	09/07/2015	06/12/2022	06/12/2022	
Init.	NV	PM	AZ	

Material				Material				
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity	
1	YP09307DP	26435...	1	32		Sterknop M6	2995334	2
2	YP00372	18117...	1	33		M12 DIN 127 A	1001572	2
3	Y004067	15072...	1	34		Handle and Line	1909113	1
4	Y004072	15073...	2	35		XP00071	2990040	2
5	YP09308DP	2643604	2	36		Screw D4x12	1990816	4
6	---	1611704	1	37		X004865	1003901	2
7	Front wheels	1916560	2	38		---	1994587	1
8	M12x60	1993419	2	39		X002335	1993411	1
9	Forks	1991447	2	40		M5x10 DIN 916/45H	1001485	2
10	Rear wheels	1914606	2	41		L=980mm	1003279	2
11	M10x80 DIN931	1001457	2			/	/	/
12	M10 DIN985	1001683	2	42		---	1996192	2
13	M6x75 DIN 912/8.8	1993413	2	43		---	---	---
14	X002306	1502304	2	44		XP00028	2990012	2
15	X002443	1517104	2	* 45		---	PP00381	1
16	X002448	1517200	2	* 46		LOGO VERMEIREN	T001002	1
17	D6,4 DIN 125 A	1001564	16	47		M6x35 DIN931	2995333	2
18	---	1003558	2	* 48		---	W620150	2
19	M6 DIN985	1001680	10	* 49		---	W620151	2
20	M6x40	1006111	2	* 50		---	1003558	2
21	M6x60 DIN 912/8.8	1005957	4	51		M8x76 DIN912	1998993	2
22	M6 DIN125A	1000973	4	52		M8	1998994	2
23	XP00099	2990067	4	*Not depicted				
24	D10,5 DIN 125 A	1001558	6					
25	M6x16	1008459	10					
26	D4x30	1001545	4					
27	X002305	1993120	2					
28	D5x25	H016007	2					
29	X005513	1992600	2					
30	XP00108	GP00604	2					
31	---	1993417	1					

EXPLODED VIEW 276				Plan No.
2/2				E000994/H
Scale:	Drafted	Verified	Confirmed	
Date	09/07/2015	06/12/2022	06/12/2022	
Init.	NV	PM	AZ	